

Red HDPE → Natural HDPE Colour Change on 6L Accumulator EISA

The Process



Poly Clear HMEX (HT)



Fourth bottle acceptable
3 rejects



Results

Normal results:

- Purge with HDPE
- 600 kg
- 4 hours

Aquapurge results:

- Poly Clear HMEX (HT)
- 60 kg
- 1 hour
- 3 rejects

Request your FREE trial today

Aquapurge Ltd, Unit 2, Argent Trade Park, Pump Lane, Hayes, Middlesex
UB3 3NB

+44 (0) 20 8813 7990 enquiries@aquapurge.com

www.aquapurge.com